

# PICO System

## Competency checklist performance criteria

Name \_\_\_\_\_ Date \_\_\_\_\_

Evaluator \_\_\_\_\_

| PICO dressing application  | Completed | Not completed |
|--|-----------|---------------|
| Discuss indications, contraindications, warnings and precautions for use of NPWT   |           |               |
| Troubleshoot loss of seal at dressing site (if applicable)   |           |               |
| Demonstrate closed incision dressing application   |           |               |
| Demonstrate open wound dressing application and use of appropriate fillers (black foam and antimicrobial gauze) when necessary |           |               |
| PICO device  | Completed | Not completed |
| Discuss an understanding of the visual alarms  |           |               |
| Identify the battery life associated with the device and therapy time  |           |               |
| Demonstrate ability to load batteries inside of pump and attach to the dressing  |           |               |
| Demonstrate turning the device on by pressing the orange button to start delivery of NPWT                                      |           |               |

### NPWT Clinical and Product Support

- Clinical Hotline 1-800-876-1261
- Reimbursement Hotline 1-888-705-0061

Additional information available at: [www.possiblewithpico.com](http://www.possiblewithpico.com)

#### Important Safety Information:

The PICO pumps contain a MAGNET. Keep the PICO pumps at least 4 inches (10 cm) away from other medical devices at all times. As with all electrical medical equipment, failure to maintain appropriate distance may disrupt the operation of nearby medical devices.